

## Introduction

This document describes the PHI collection and storage of Cognivue Clarity and Thrive. It also provides details on the security measures on the device.

## PHI Collection and Storage

Cognivue Clarity and Thrive collect and store the following information on a per-test basis:

- Patient's first name, middle initial, and last name.
- Patient's date of birth
- Patient's gender
- Patient's subtest scores and domain scores.
- (Optionally) Patient's email.

Cognivue Clarity and Thrive store information on the device in a Postgres database that is only accessible from on-device processes.

As part of the screening process, the device also generates a report which contains this information.

Administrators can set up a network file share to access the directory where the reports are stored.

A device administrator can also configure the device to create an HL7 message that holds the listed PHI.

## Device Security

As part of our responsibility to comply with the FDA's guidance on medical device cybersecurity and HIPAA regulations, the device has several security features in place.

- The device is managed and monitored by the AWS [Systems Manager](#) tool. This tool enables our team to resolve issues remotely; Systems Manager uses HTTPS to communicate.
- The on-device firewall blocks all incoming and outgoing connections except for the ones needed for the file share. (Administrators can find details at [Firewalling Samba](#))
- Cognivue Clarity/Thrive uses HTTPS for all outgoing and incoming network traffic.
- Access to the underlying OS (Ubuntu 18.04 LTS) is not permitted except for limited operations.
  - These include network and printer configuration.
  - The device can use a printer connected by USB or a network at will.
- The OS receives security updates automatically.
- The system has an AV solution, [ClamAV](#), installed and running.

- The folder share requires a password that the device administrator (Admin) can change.
- Device administrators are the only user that view reports on the device.
- Operator users can view the test reports immediately after the test concludes but cannot view older reports on the device.
- Device administrators can enable user account password expiration.
- The device can hide all reports from tests administered before a specific date.
  - This action does not delete the database results; however, it removes the reports from the Reports folder and hides the corresponding test session from the UI.