



# User Manual

# Cognivue Thrive® User Manual

## 1 Contents

2	About This Guide.....	4
3	About Cognivue Thrive® .....	4
3.1	Intended Purpose.....	4
3.2	Indications for Use .....	4
3.3	Clinical Reference Information .....	4
4	Device Description .....	5
5	Power Management .....	6
5.1	Indicator Lights.....	6
5.2	Screen Sleep.....	6
5.3	Power Alert .....	6
6	Recommendations for Use .....	7
6.1	General Testing Location .....	7
6.2	Cleaning and Disinfecting Cognivue Thrive® .....	7
6.3	Use of Personal Protective Equipment .....	8
7	Network Connectivity .....	8
7.1	Basic Wired Connection .....	8
7.2	Basic Wi-Fi Configuration .....	9
7.3	Connection Details.....	9
8	Security .....	10
8.1	Login.....	10
8.2	User Management .....	11
8.3	User Creation .....	12
8.4	User Modification .....	12
8.5	Password Reset .....	13
8.6	Auto-Logout .....	13
8.7	Password Expiration.....	14
8.8	Disabling Accounts .....	14
9	System Settings.....	15
9.1	Client Code .....	15
9.2	Video Settings .....	16
10	Report Settings.....	17
10.1	Report Emailing.....	17

# Cognivue Thrive® User Manual

10.2	Report Printing .....	18
10.3	Automatic Printing .....	18
10.4	Report Composition .....	19
10.5	Practice Information .....	19
11	Portal Menu .....	19
12	Patient Management .....	20
12.1	Pre Entering Patients .....	20
13	Assessments .....	20
13.1	Run Test (Main Screen) .....	21
13.2	Run Test (Patient List Screen) .....	22
13.3	Video Volume Control .....	23
14	Assessment Reports .....	24
14.1	Practitioner Report .....	24
14.2	Patient Report .....	25
14.3	Cogniwell Letter .....	26
15	Accessing Reports .....	27
15.1	Viewing Reports .....	27
15.2	Reviewing Reports .....	27
15.3	Printing .....	27
15.4	Emailing .....	28
15.5	Drive Mapping .....	29
15.5.1	Setting the SMB Share Password .....	29
15.5.2	Reports Mapping .....	29
15.5.3	CSV Mapping .....	30
16	Technical Specifications .....	31
17	Safety Precautions .....	31
18	Product Safety Standards .....	32
19	Product Safety Symbols .....	33

# Cognivue Thrive® User Manual

## 2 About This Guide

This User Manual contains information that describes the operation of Cognivue Thrive®. Additional information can be found at [www.cognivue.com](http://www.cognivue.com). Cognivue Support personnel can be reached via at 1-585-433-2992 from 9-5 pm ET or by emailing [support@cognivue.com](mailto:support@cognivue.com).

UM-420-U Cognivue Thrive® User Manual 3/22/2023



© 2023 Cognivue, Inc. | 7911 Rae Blvd | Victor, NY 14564 | 585.203.1969

## 3 About Cognivue Thrive®

Cognivue Thrive® is a computerized cognitive assessment device to evaluate cognitive function objectively, quantitatively, and reliably. Cognivue Thrive® relies on a series of 6 subtests. Using the results from these subtests, Cognivue Thrive® creates a short, easy-to-interpret report. The report shows three cognitive domain scores: Memory, Visuospatial, and Executive Function. Below these are two metrics on Reaction Time and Speed Processing. Normative ranges appear below the visual of each cognitive domain score and metric. These segment into classifications of good, moderate, or poor ability. These cognitive scores and classifications provide a baseline for future comparison and enable discussion between the clinician and the test-taker about cognitive health.

### 3.1 Intended Purpose

Cognivue Thrive® is indicated for use as an adjunctive tool for evaluation cognitive function in ages 18 and over.

### 3.2 Indications for Use

Cognivue Thrive® is indicated for use as an adjunctive tool for evaluating cognitive function in ages 18 and over. It is not a stand-alone diagnostic tool. The device results are to be assessed and interpreted by a licensed clinician.

### Warnings

- 1) Cognivue Thrive® does not identify the presence or absence of clinical diagnoses.
- 2) Cognivue Thrive® cannot be used as a stand-alone diagnostic.
- 3) Cognivue Thrive® is intended to be used in clinical professional offices or similar settings.

### 3.3 Clinical Reference Information

Cognivue Thrive® is based on the same technology as Cognivue® Advanced and Cognivue® Clarity. Cognivue Thrive uses a subset of the Cognivue Advanced® algorithms to provide a similar breakdown of cognitive scores in a shorter duration test.

Cognivue Thrive® algorithms have been compared vs. Cognivue Advanced® in a clinical trial with 100 subjects. The correlation between Advanced and Thrive was in the same category of Advanced vs.

# Cognivue Thrive® User Manual

Advanced, which was statistically significant (moderate positive 0.57), and Thrive's test reliability was high positive (0.74).

Additional information is available in the following references:

- Cognivue Advanced De Novo Clearance, DEN130033
- Diego Cahn-Hidalgo, Paul W Estes, Reina Benabou, "Validity, reliability, and psychometric properties of a computerized cognitive assessment test (Cognivue®)," World Journal of Psychiatry 2020 January 19; 10(1): 1-11.
- Cognivue Clinical Analysis Report CAR-401-A.

## 4 Device Description

The **Cognivue Thrive®** is a portable, battery-powered, laptop-like device that includes a wireless keyboard with a touchpad, a rotary CogniWheel®, and a CogniCover®. The CogniWheel® is the sole patient input device. The CogniCover® is designed for privacy and to reduce distraction and light reflections on the computer screen. The CogniCover® and CogniWheel® are designed to fold for portability.




Figure 1 Cognivue Device (Front)

## 5 Power Management

The Cognivue Thrive®, much like a laptop, is designed to be battery-powered during regular use but can also be used when the provided charging adapter is plugged in.

The battery sustains 4 hours of intensive, off charger usage but can last for over 8 hours typical usage. The battery recharges in approximately 4 hours if depleted completely.

### 5.1 Indicator Lights

There are two indicator lights on the right side of the device.

	<i>Left Light</i>	<i>Right Light</i>
<i>Device is unpowered</i>	Off	Off
<i>Device is discharging</i>	Orange	Off
<i>Device is charging</i>	Green/Off	Orange
<i>Plugged in and fully charged</i>	Green/Off	Green

Table 1 Charging Lights



Figure 2 Cognivue Device (Rear)

USB (for Keyboard)

Ethernet

Battery Charger

### 5.2 Screen Sleep

When the Cognivue Thrive device is not in use for 10 minutes, the screen will go into sleep mode. When the screen is asleep, it can be awakened by clicking or pressing a key on the keyboard.

### 5.3 Power Alert

When the device is low on battery charge and not plugged in, if an attempt is made to start a testing session, e.g., upon clicking the **Run a Test** button on the main screen, it will display an alert and not allow the testing session to begin. Plugging in the device will allow tests to be administered. In order to avoid running out of power during a test, please charge the device for at least one hour before again administering tests on battery power.

# Cognivue Thrive® User Manual

## 6 Recommendations for Use

### 6.1 General Testing Location

When possible, choose a test location that reduces distractions (e.g. an exam room, office, or conference room rather than a busy hallway).

Put the device on a flat, stable surface that is ergonomically comfortable for the patient. The patient should be able to comfortably reach and manipulate the CogniWheel® and view the computer screen. Avoid placing the Cognivue Thrive® in a location where there may be distractions behind the device which are visible to the patient.



Figure 3 Example Setup

**DO NOT place the device on the patient's lap for the test.** This does not provide a stable surface during the administration of the test. This interferes with a patient's ability to take the test.

### 6.2 Cleaning and Disinfecting Cognivue Thrive®

Cognivue recommends sanitizing the Cognivue Thrive® before and after administering a test using one of the following procedures:

- **Preferred:** Use sanitizing or disinfectant wipes to thoroughly wipe the CogniWheel® and the base of the device where the patient palm rests. **DO NOT** wipe the screen!
- Please turn the power off on the Cognivue Keyboard using the appropriate slider switch prior to cleaning/transport; this prevents unintentional keypress combinations on the device.
- In the event wipes are not available, spray a sanitizer onto a paper towel, being careful not to oversaturate the paper towel, and thoroughly wipe the CogniWheel® and the base of the device where the patient palm rests. **DO NOT** wipe the screen, and do NOT spray sanitizer or disinfectant directly onto the device!

# Cognivue Thrive® User Manual


## 6.3 Use of Personal Protective Equipment

In addition to sanitizing the device, Cognivue recommends that operators and patients wear a disposable mask while interacting with the device. If it is not possible to sanitize the device between patients, Cognivue recommends operators and patients use disposable gloves.

## 7 Network Connectivity

Cognivue Thrive® requires connection to a qualified network for initial setup and configuration but thereafter can be used to administer assessments offline. However certain features like Portal Sync, Printing, and Drive Mapping will be unavailable until the device is reconnected to a qualified network. In addition to these features, continued connection to qualified networks allows Cognivue to deploy both general purpose and security patches to the devices.


### 7.1 Basic Wired Connection

1. Connect one end of a standard network cable to an active network jack at your location
2. Connect the other end to the network port on the back of the Cognivue Thrive® device
3. If not automatically connected, follow the instructions above, selecting **Wired Connection 1**.
4. Success is indicated when the icon in the top left corner is the  Network icon.



# Cognivue Thrive® User Manual

## 7.2 Basic Wi-Fi Configuration

1. Click on the  Network icon in the status bar at the upper left corner of the screen
2. From the drop-down menu, select the available Wi-Fi network to which you wish to connect
  - a. If there are a large number of networks available, you may see a **More networks** option which you can check if you do not see your network listed.




Figure 4 Network Menu

3. Enter the Wi-Fi password for the selected network
4. Click **Connect**



Figure 5 Authentication Prompt

automatically connect when possible.



Figure 6 Network Selection Menu

## 7.3 Connection Details

When logged on to a Cognivue Thrive® device, administrators can view and test the current network configuration information on the Network Configuration window

To open the **Network Configuration** window

# Cognivue Thrive® User Manual

1. Log in
2. Select **Device Administration** >> **Settings**
3. Click the **Network** tab

The **Network Configuration** window contains network configuration properties and their values

- The **Refresh** button refreshes the displayed connection values
- The **Test Connection** button confirms that the device has network connectivity
- The **Advanced Configuration** button opens settings for more advanced networking options

Name	Signal	Security
ATT-WIFI-0553	99	WPA2
CogniGuest	95	WPA2

Network: CogniNet

IP Address: 192.168.0.122

Netmask: 255.255.255.0

Gateway: 192.168.0.1

MAC Address: DC-1B-A1-4E-88-5B

Password:

Connect

Refresh

Test Connection

Advanced Configuration

Figure 7 Network Configuration Screen

Advanced Configuration is intended for use by network administrators and/or Cognivue Support Staff.

The most common usage of this functionality is outlined in the guide for Static IP Configuration

## 8 Security

### 8.1 Login

The user **Login** window is how you log in to a device.

It is how you log in to the device. Logging in with an individual account ensures an audit trail for viewing patient reports and running assessments.

Username:

Password:

Login

Shutdown

Figure 8 Login Screen

The **Shutdown** button appears on the **Login** window can be used by anyone to power down the device. Clicking the **Shutdown** button will prompt confirmation before shutting down the device.

**Note:** This is the preferred way to shut down the device. The physical power button should be used only if necessary.

# Cognivue Thrive® User Manual

## 8.2 User Management

Each Cognivue Thrive® user should have their own user account to ensure audit trail and safety of patient data.

User accounts can be assigned as an Administrator or an Operator on the **User Management Screen**.

	Administrator	Operator
Setting up new accounts	X	
Delete accounts	X	
Changing accounts	X	
Configure network connection	X	X
Entering patient data & running a test	X	X
Running a test with an intro video	X	X
Viewing past test reports	X	

Table 8-1

To open the **User Management Screen**

1. Log in
2. Select **Device Administration >> Settings**
3. Click the **Users** tab





Figure 9 User Management Screen

# Cognivue Thrive® User Manual

## 8.3 User Creation

From the **User Management Screen**

1. Click **Create**
2. Fill out form
3. Click **Save**
4. Click **Yes** to confirm when prompted




The 'User Creation' form is a light gray box with a title bar containing a user icon and the text 'User Creation' in blue. A red 'X' icon is in the top right corner. Below the title bar, a subtitle reads: 'Complete the form with the information for the new user and click Save.' The form contains several labeled input fields, each preceded by a red asterisk indicating it is required: 'Username:' with an empty text box; 'Password:' with a text box containing 'password'; 'Confirm Password:' with a text box containing 'confirm password'; 'First Name:' with a text box containing 'first name'; 'Last Name:' with a text box containing 'last name'; 'Email:' with a text box containing 'email'; 'Confirm Email:' with a text box containing 'confirm email'; and 'Role:' with a dropdown menu showing 'OPERATOR'. At the bottom left, a red asterisk is followed by the text '\* Required field'. At the bottom right, there is a 'Save' button.

Figure 10 User Creation Screen

## 8.4 User Modification

From the **User Management Screen**

1. Select a user from the list
2. Click **Edit**
3. Make the desired changes
4. Click **Save**
5. Click **Yes** to confirm when prompted



The 'User Modification' form is a light gray box with a title bar containing a user icon and the text 'User Modification' in blue. A red 'X' icon is in the top right corner. Below the title bar, a subtitle reads: 'Make any desired modifications to the information for the selected user and click Save. To set the user password, click the Set Password button.' The form contains several labeled input fields, each preceded by a red asterisk indicating it is required: 'Username:' with a text box containing 'ADMIN1'; 'First Name:' with a text box containing 'Sam'; 'Last Name:' with a text box containing 'Smith'; 'Email:' with a text box containing 'ssmithe@company.co'; and 'Role:' with a dropdown menu showing 'ADMINISTRATOR'. At the bottom left, there is a 'Set Password' button. At the bottom right, there is a 'Save' button. At the bottom left, a red asterisk is followed by the text '\* Required field'.

Figure 11 User Modification Screen

## 8.5 Password Reset

From the **User Management Screen**

1. Select a user from the list
2. Click **Set Password**
3. Fill out the Set Password form
4. Click **Save**
5. Click **Yes** to confirm when prompted


A screenshot of a web form titled "Set User Password" with a user icon and a close button (X) in the top right corner. The form contains the instruction "Please enter your new password". There are two input fields: the first is labeled "\*Password:" and is empty; the second is labeled "\*Confirm Password:" and contains the placeholder text "confirm password". A red asterisk label "\* Required field" is positioned to the left of the "Save" button at the bottom right of the form.

Figure 12 Password Reset

## 8.6 Auto-Logout

The device will be logged out after 30 minutes of inactivity to prevent unauthorized access.

# Cognivue Thrive® User Manual

## 8.7 Password Expiration

Password expiration forces users to change their account passwords after a set number of days. Password expiration applies to all user accounts but is disabled by default.

Users are warned about password expiration for three days before it expires. Once expired, the user will be forced to change their passwords before they can log in.




Figure 14 Password Expiration Section




Figure 14 Password Expiration Prompt

To enable/disable password expiration:

1. Navigate to **Device Administration >> Settings**
2. Click the **Enable Password Expiration** checkbox
3. When this feature is enabled, a password lifetime must be specified
4. Enter the number of days in the text field labeled **Password Expires after**

## 8.8 Disabling Accounts

To maintain an audit trail, User Accounts cannot be fully removed from the system. Disabling prevents the user from logging in. Disabled accounts can be reactivated in the future.

From the **User Management Screen**

1. Click on a user in the list
2. Click the **Disable User** button
3. Click **Yes** to confirm when prompted




Figure 15 Disable User Prompt

## 9 System Settings

The Cognivue Thrive® device includes several configurable options. The ability to make these configuration changes is available only to Administrator users on the **System Settings** screen

The System Settings Screen is on the System tab of the Device Administration page

To open the **System Settings** Screen

1. Log in
2. Select **Device Administration**  
    >> **Settings**
3. Click the **System** tab




Figure 16 System Settings Screen

### 9.1 Client Code

The Cognivue Thrive® device provides the option to include a custom code on all tests for purposes of categorization. Any sequence of up to 100 characters can be entered for each patient test session.

A user account with the Administrator role can enable it by checking the **Enable Client Code** checkbox on the **System Settings** screen. Making this selection will add the Client Code text box to the **Assessment Setup** screen before each assessment, where operators can enter the client code.




Figure 18 System Settings Screen (Client Code)




Figure 17 Assessment Setup Screen (Client Code)

## 9.2 Video Settings

This screen allows Administrators to ensure certain videos will always be shown to patients. If the option is unchecked here, the choice is left to the Operator who sets up the test for the patient.

- Always Play Intro Video → Play the introduction video before every test
- Always Play Instructional Videos → Plays short instructional videos during the test
- Always Play Post Test Video → Plays a congratulatory video when the test completes




Figure 20 Default Video Setting Screen

Figure 20 Assessment Setup Screen (Videos)



## 10 Report Settings

The Thrive reports are configurable to meet the needs of the practice using the device. These configuration changes are available only to Administrator user accounts on the **Report Configuration** screen.

### 10.1 Report Emailing

Cognivue Thrive® provides the ability to send the test report via an encrypted email to a patient following completion of a test. This feature is enabled by default. An Administrator can disable this functionality by unchecking the **Enable Report Emailing** checkbox.




Figure 22 Report Settings Screen




Figure 23 Assessment Setup Screen (Email)

When enabled, the option for the patient to receive a copy of the test report via encrypted email will be provided on the **Assessment Setup** Screen.




Figure 21 Default Report Settings Screen

# Cognivue Thrive® User Manual

## 10.2 Report Printing

Cognivue Thrive® provides the ability to print the test report to a configured, connected printer. Printing test reports from Cognivue Thrive® is enabled by default. An Administrator can disable this functionality by unchecking the **Enable Report Printing** checkbox.




Figure 25 Report Settings Screen




Figure 24 Report Viewer Screen

When enabled, a **Print** button will be displayed on the top navigation bar of the **Report Viewer** screen.

## 10.3 Automatic Printing

Cognivue Thrive® provides the ability to automatically print the associated report upon completion of a test. Automatic printing of test reports is disabled by default. An Administrator can enable this functionality by checking the **Auto Print** checkbox




Figure 27 Reports Settings Screen

# Cognivue Thrive® User Manual

## 10.4 Report Composition

To open the **Report Settings Screen**

1. Log in
2. Select **Device Administration >> Report Settings**
3. Click the **Composition** tab



Figure 28 Report Composition Screen

Report Composition allows the inclusion of the Cogniwell Information Sheet and picking one or both reports.

**WARNING** Only the report generated will be available.

## 10.5 Practice Information

To open the **Practice Information Screen**

1. Log in
2. Select **Device Administration >> Report Settings**
3. Click the **Practice** tab




Figure 29 Practice Info Screen


## 11 Portal Menu

For details about how to setup and utilize the Portal see the **Portal User Manual**

## 12 Patient Management

To open the **Patient List** Screen

1. Log in
2. Click the **Patient** button



The screenshot shows the 'Patient List' screen with a title bar, a close button (X), and two buttons: '+ Create' and '▶ Start Test'. Below these is a table with three columns: 'Patient Name', 'Date of Birth', and 'Last Test Taken'. The table contains three rows of data, with the third row highlighted in blue.


Patient Name	Date of Birth	Last Test Taken
Jessica L Anque	02/11/1955	08/24/2022 2:14 PM
John D Doe	12/13/1988	08/24/2022 12:53 PM
Dave X Mage	12/22/2000	

Figure 30 Patient List Screen

### 12.1 Pre Entering Patients

Patient demographic information can be pre-filled to speed up the test setup process

1. Click on the **Create** button to display the **Patient Creation** Screen
2. Fill out form
3. Click **Save**
4. Click **OK** to confirm when prompted



The screenshot shows the 'Patient Creation' screen with a title bar, a close button (X), and a user icon. Below the title bar is a message: 'Complete the form with the information for the new patient and click Save.' The form has several fields: 'First Name\*' (John), 'MI\*' (X), 'Last Name\*' (Smith), 'Date of Birth\*' (01 / 02 / 2003), 'Sex\*' (Male), and 'Email Address'. A 'Save' button is at the bottom right. A red asterisk indicates required fields.

Complete the form with the information for the new patient and click Save.

First Name\* MI\* Last Name\*

John X Smith

Date of Birth\*

01 / 02 / 2003

Sex\* Male

Email Address

\* Required field

Save

Figure 31 Patient Creation Screen

## 13 Assessments

Test administration is the core feature of Cognivue Thrive®. There are two paths to starting an assessment: From the **Main** Screen and from the **Patient List** Screen.

# Cognivue Thrive® User Manual

## 13.1 Run Test (Main Screen)

1. From the **Main Screen** Click the **Run a Test** button




Figure 13-1 Main Screen (Administrator)

2. Fill out patient data (Name, \*sex, date of birth)
3. Select the patient's preferred language for the test (English or Spanish)
4. Enter any additional information as desired
5. Select appropriate video options or use the defaults
6. Click **Run Test**

The image shows the 'Patient Info' form. It contains fields for 'First Name\*' (John), 'MI\*' (X), and 'Last Name\*' (Smith). Below these is the 'Date of Birth\*' field with a date picker set to 01 / 01 / 2001. There is a 'Sex\*' dropdown menu set to 'Male' and a 'Test Language' dropdown menu set to 'English'. At the bottom, there are expandable sections for 'Email' and 'Video Options'. At the very bottom are three buttons: 'Run Test', 'Cancel', and 'Test Audio'. Blue arrows point from the steps 2 through 6 to the corresponding fields and buttons in the form.

Figure 13-2 Assessment Setup Screen

7. The patient should turn the CogniWheel when ready to begin





Figure 32 Start of the Cognivue Assessment

## 13.2 Run Test (Patient List Screen)

From the **Patient List** Screen

1. Select a patient
2. Click the **Start Test** button




The Patient List screen shows a table with three columns: Patient Name, Date of Birth, and Last Test Taken. There are three patients listed: Jessica L Anque, John D Doe, and Dave X Mage. The 'Start Test' button is highlighted in orange.

Patient Name	Date of Birth	Last Test Taken
Jessica L Anque	02/11/1955	08/24/2022 2:14 PM
John D Doe	12/13/1988	08/24/2022 12:53 PM
Dave X Mage	12/22/2000	

Figure 35 Patient List Screen

3. Fill out patient data (Name, \*sex, date of birth)
4. Select the patient's preferred language for the test (English or Spanish)
5. Select appropriate video options or use the defaults
6. Enter any additional information as desired
7. Click **Run Test**



The Assessment Setup screen contains fields for Patient Info. It includes input fields for First Name, MI, and Last Name, a Date of Birth field, a Sex dropdown menu, and a Test Language dropdown menu. There are also expandable sections for Email and Video Options, and buttons for Run Test, Cancel, and Test Audio.

Figure 35 Assessment Setup Screen

8. The patient should turn the CogniWheel when ready to begin




Figure 35 Start Assessment Screen

## 13.3 Video Volume Control

Before starting Cognitive Assessments test the audio using the **Test Audio** button on the **Assessment Setup** Screen




Figure 36 Assessment Setup Screen

In order to control the volume of video playback, the volume control keys on the keyboard can be used.



Table 13-1 Volume Keys

**Note:** The **Mute** button sets the volume to its lowest setting. To unmute, it is necessary to use the volume-up key to increase the volume to the desired level.

# Cognivue Thrive® User Manual

## 14 Assessment Reports

### 14.1 Practitioner Report

This document is intended for the practitioner. It contains the test scores and normalized ranges. Additional information is provided to the practitioner to help with explaining scores to the patient.




Figure 37 Practitioner Report



# Cognivue Thrive® User Manual

## 14.2 Patient Report

This document is intended for the patient. It contains the test scores and normalized ranges. Additional information is provided to define the domains in patient-friendly language.




Figure 14-2 Patient Report

# Cognivue Thrive® User Manual

## 14.3 Cogniwell Letter

It provides information about Cognivue's Cogniwell program



**Take control of your memory and brain health**

You have the power to optimize your cognition by adopting a combination of key healthy lifestyle habits and being aware of your overall wellness. Studies show that taking proactive steps may help reduce your risk of memory loss and other forms of cognitive decline. Here are evidence-based recommendations to proactively manage lifestyle-related risk factors and certain medical conditions that may affect cognition. For more information, please contact one of our Cogniwell Coaches at [wellness@cognivue.com](mailto:wellness@cognivue.com) or (844) 654-1009



**Stay Mentally Fit**

**Exercising your mind can help your brain perform at its best**

Participating in cognitively stimulating activities may help sharpen certain thinking skills that tend to diminish with age. Have fun engaging your brain by reading books, playing games, and doing puzzles.



**Maintain Connections**

**Being socially active is brain friendly**

Studies have shown that increased social activity is linked to a lower rate of cognitive decline. To stay social, connect with friends and family, or consider volunteering for a local organization.



**Exercise Regularly**

**Exercising your body can help keep your brain in shape**

Research supports a strong link between physical activity and brain health. It is recommended to get at least 150 minutes of moderate aerobic activity a week. Any activity that gets your heart pumping counts as aerobic activity, such as walking, dancing, and biking.



**Eat Healthy**

**Eating right feeds a healthy body and mind**

Food that is good for your overall health is good for your brain. Talk to your physician about whether you are consuming a balanced diet that takes into consideration appropriate numbers of calories and nutrients.



**Sleep & Relax**

**Getting a good night's sleep rests your body and mind**

Insufficient sleep can not only make a person irritable, it can affect memory and decision making. Healthy adults need between 7 and 9 hours of sleep per night.



**Healthy Hearing**

**Hearing health is important to your cognitive health**

Hearing impairment has been recognized as a modifiable contributor to cognitive decline when identified in mid-life. See your hearing healthcare professional to help maintain your hearing health.



**Optimal Vision**

**Manage or prevent vision impairment**

Research has found the relationship between cognitive decline and vision impairment may be a modifiable risk factor. To manage or prevent vision impairment, schedule regular visits with your vision healthcare provider.



**Medication Management**

**Medication management is important to your overall health**

Taking your medicine as prescribed is important for overall health. Reviewing your medications with your healthcare provider or pharmacist is an integral part of medication management.

For a free Consultation, contact our Wellness Coach at (844) 654-1009

 **cognivue**  
cognivue.com

Disclaimer: Cognivue Clarity/Thrive® is indicated for use as an adjunctive tool for evaluating cognitive function. It is not a stand-alone diagnostic tool and does not identify the presence or absence of clinical diagnoses. The device results are to be assessed and interpreted by a licensed clinician. Cognivue, Cognivue Thrive and Cognivue Clarity are trademarks or registered trademarks of Cognivue in the US and/or other countries. © Cognivue. All rights reserved. CG-1003 (07/2022)

Figure 14-3 Cogniwell Letter


## 15 Accessing Reports

### 15.1 Viewing Reports

After a Cognitive Assessment is completed When the test has been completed the report can be reviewed by going to the **Report List** Screen.

### 15.2 Reviewing Reports

Administrators can review the test report by Selecting the **Reports** button from the Main Menu Screen. A list of available reports will be presented. After selecting the desired report from the list, click the **View Report** button to display the report.



The screenshot shows the 'Reports List' screen with a 'View Report' button and a table of reports.

Last Name	First Name	Middle Initial	Birthdate	Test Time	Language	Sync
Smith	John	X	01/01/2001	09/14/2022 10:46 AM	ENGLISH	Y
Smith	John	X	01/01/2001	09/13/2022 2:29 PM	ENGLISH	Y
Smith	John	X	01/01/2001	09/13/2022 2:04 PM	ENGLISH	Y

Figure 38 Report List Screen

### 15.3 Printing

Test reports may be printed from the Thrive device whenever a report is opened in the report viewer. The print feature prints to the device's default printer.

The Default printer can be set from any screen

1. Click the printer icon using the menu bar in the upper left corner of the screen
2. Right click on the preferred Printer
3. Click **Set As Default**

A green checkbox indicates that the printer has been set as default




Figure 40 Printer Management Screen

## 15.4 Emailing

Emailing reports from the device requires consent on a per assessment basis. The checkbox on the start assessment screen or the re-entry of a pre-existing patient's email serves as this consent.

Patient Info

▼ Email

☐ Email Report

By checking this box, you are agreeing to receive emails from Cognivue.

▶ Video Options

Run Test Cancel Test Audio

Figure 42 Email Report Checkbox

Patient Info

▼ Email

☒ Email Report

By checking this box, you are agreeing to receive emails from Cognivue.

Email Address\*

jsmith@cognivue.com

Confirm Email Address\*

Email addresses must match

▶ Video Options

Run Test Cancel Test Audio

Figure 42 Email Report Confirmation

If the email cannot be sent after the assessment due to network interruption the email can be resent once the network is restored via the **Email Menu** on the **Main Screen** for an administrator user.




Figure 43 Email Report Menu

# Cognivue Thrive® User Manual

## 15.5 Drive Mapping

### 15.5.1 Setting the SMB Share Password

The SMB password can be set and reset on device by an administrator user.

From the **System Settings** Screen

1. Click the **SMB Share** Button




Figure 46 System Settings Screen

2. Enter a password that meets all the requirements

The screenshot shows the 'Change SMB Share Password' dialog. It has two password input fields: 'New Password:' and 'Confirm Password:'. Below the 'New Password' field, there are four requirements listed in green: 'At least 8 characters long', 'At least 1 lowercase letter', 'At least 1 capital letter', and 'At least 1 number'. A 'Change Password' button is at the bottom right. A blue arrow points from the 'New Password' field to the next step.

Figure 46 Change SMB Password

3. Click **Change Password** Button

The screenshot shows the 'SMB Change Confirmation' dialog. It has a title bar 'Change SMB Password' and a 'Confirmation' section with a question mark icon. The text asks 'Are you sure you want to change the SMB share password?'. There are 'Yes' and 'No' buttons at the bottom. A blue arrow points from the 'Yes' button to the next step.

Figure 46 SMB Change Confirmation

4. Click to confirm the change

### 15.5.2 Reports Mapping

This mapping provides access to the report pdfs sorted into folders based on the name and date of birth of the patient.



# Cognivue Thrive® User Manual

To map this drive, follow the appropriate Drive Mapping Instructions available at [www.cognivue.com/support](http://www.cognivue.com/support).

The information you will need for this share is

- Ip Address: Retrieved from the Network tab
- Mapping folder: Reports (e.g., \\192.168.0.48\Reports)
- Username: cogaccess
- Password: Set by device administrator

## 15.5.3 CSV Mapping

The CSV store provides access to csv files with the condensed data from all patient reports. These exports are generated by clicking the **Data Export** menu and selecting the **CSV Export** menu item.




Table 15-1 Main Screen (CSV Export)

The data included is:

- Last Name
- First Name
- Middle Initial
- Patient ID
- Sex
- Date of Birth
- Test Start Date-Time
- Serial Number
- Language
- Intro Video
- Instructional Videos
- Post Test Video
- Visuospatial
- Executive Function
- Memory
- Reaction Time (ms)
- Processing Speed (ms)

# Cognivue Thrive® User Manual

To map this drive, follow the appropriate Drive Mapping Instructions available at [www.cognivue.com/support](http://www.cognivue.com/support).

The information you will need for this share is

- Ip Address: Retrieved from the Network tab
- Mapping folder: Csv (e.g., \\192.168.0.48\Csv)
- Username: cogaccess
- Password: Set by device administrator

## 16 Technical Specifications

<b>Operating Environment</b>	10 to 30 °C (50 to 86 °F) 10% to 85% Relative humidity, non-condensing 70 to 106 kPa Atmospheric Pressure
<b>Transport and Storage Environment</b>	-20 to 50 °C (-4 to 122 °F) 10 to 85% Relative humidity, non-condensing 50 to 106 kPa Atmospheric Pressure
<b>Power Input</b>	19V DC, 3.5 A Nominal
<b>Weight</b>	4.2 kg (9lbs 2oz)
<b>Wireless Communications</b>	Intel Dual Band Wireless-AC 3168 20 dBm max output power for 2.4 GHz 802.11b/g/n 23 dBm max output power for 5.0 GHz 802.11a/n/ac

## 17 Safety Precautions

Observe the following precautions to ensure safety and prevent damage to Cognivue Thrive®:

- Do not place on uneven or unstable surfaces.
- Do not make any modifications or alterations to the device.
- Do not connect unauthorized peripherals to the device.
- Do not open the device cover or modify the device. There are no customer serviceable parts.
- Do not apply weight onto the device.
- Do not connect a multiple socket outlet or extension cords to Cognivue Thrive®.
- Do not install the unit such that the power cord is a tripping hazard.
- Do not expose to liquids, rain, or moisture except in routine cleaning of surfaces.
- This product contains rechargeable batteries. At end of use, return the device to Cognivue or dispose according to local guidelines.
- Only use the provided AC Power Adapter with this device.
- Removal of the power cord from the inlet is the mains disconnect.

**WARNING:** To avoid the risk of electric shock, this equipment must only be connected to a supply mains

# Cognivue Thrive® User Manual

with protective earth.

**WARNING:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

**WARNING:** Use of accessories, transducers, and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

**WARNING:** Portable RF communications equipment should be used no closer than 30 cm (12 inches) to any part of Cognivue Thrive®. Otherwise, degradation of the performance of this equipment could result.

**NOTE:** The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

## 18 Product Safety Standards

This device meets the following product safety standards:

**USA** ANSI/AAMI ES60601-1:2005+A2 (R2012) +A1 - Medical Electrical Equipment—Part 1: General requirements for safety and essential performance

**Canada** CAN/CSA C22.2 No. 60601-1:2014 - Medical Electrical Equipment—Part 1: General requirements for safety and essential performance (includes Amendment 1)








**EMC** IEC 60601-1-2:2014 includes EMC requirements and tests. Medical Electrical Equipment including CISPR 11:2009 + A1:2010, Group 1, Class A as shown in the table below:

EMC Emission	Parameter	Standard	Test Level/Note
	Conducted Emissions	5011:2009, CISPR 11:2009 + A1:2010	Class A
	Radiated Emissions	N 55011:2009, CISPR 11:2009 + A1:2010	Class A
	Harmonic Current	EN/IEC 61000-3-2	Class A
	Voltage Flicker	EN/IEC 61000-3	-
EMC Immunity	Electrostatic Discharge	IEC/EN 61000-4-2	15 kV air, 8kV contact
	Radiated Immunity	IEC/EN 61000-4-3	3V/m (80Mhz-2.7GHz)
	Electrical Fast Transients	IEC/EN 61000-4-4	2kV
	Surge Line to Line	IEC/EN 61000-4-5	1kV Line to Line 2kV Line to Ground
	Conducted Immunity	IEC/EN 61000-4-6	3Vrms and 6Vrms
	Magnetic Immunity	IEC/EN 61000-4-8	30 A/m
	Voltage Dips & Interruptions	IEC/EN 61000-4-11	








## 19 Product Safety Symbols

The following symbols may be used for marking on this equipment:

	Follow the operating instructions
	Type B applied part
	Serial Number
	DC input
	Manufacturer Date
	Manufacturer's name and address
	NRTL Mark (Safety Certification Mark)

The following symbols may be used for Environmental Conditions for Transport

	Atmospheric Pressure Limitation
	Temperature Limitation
	Relative Humidity Limitation
	Keep away from rain
	Handle with Care

**Manufactured by:**

Cognivue, Inc.  
7911 Rae Blvd  
Victor, NY 14564  
USA

**United Kingdom Responsible Person**

Donawa Lifescience Limited  
Aviation Business Park  
Christchurch  
BH23 6NX  
UK