

Cognivue® Portal

User Manual

Contents

- 1.0 About the Cognivue Portal3
- 2.0 Portal Description.....3
- 3.0 Recommendations for Use4
- 4.0 Getting Started & User Management.....5
 - 4.1 Configuring a Cognivue Device to Sync to the Portal5
 - 4.2 Setting up the First Portal Account.....7
 - 4.3 Granting Device Users to Access the Portal (Provisioning Users).....8
 - 4.4 Adding a New User that is not a Device User to the Portal.....9
 - 4.5 Editing a user 10
- 5.0 Portal Functionality..... 11
 - 5.1 Searching for Patients & The Patient Dashboard..... 11
 - 5.2 The Reports Page..... 13
 - 5.3 Downloading Test Data 14
 - 5.4 History of Shared Reports 15
- 6.0 Devices Page 16
 - 6.1 Accessing the Devices Page 16
 - 6.1.1 Adding a Nickname to a Device 16
- 7.0 Account Management..... 18
 - 7.1 Accessing the Account Management Page 18
 - 7.2 Profile Page..... 19
 - 7.3 Password Page 19
 - 7.4 Personal Data Page..... 20
 - 7.5 Two-Factor Authentication..... 20
 - 7.5.1 Setup Authenticator App 21
 - 7.5.2 Setup SMS (texting) 21
 - 7.5.3 Recovery Codes..... 22
 - 7.5.4 Resetting Authenticator App..... 23
 - 7.5.5 Disabling Two-Factor Authentication 25
- 8.0 Portal Support..... 28
 - 8.1 Accessing the Portal Support Page..... 28
 - 8.2 Submit a Request..... 30
 - 8.2.1 Support Request 30
 - 8.2.2 Enhancement Request..... 31

About This Guide

This User Manual contains information that describes the access and operation of the Cognivue Portal. Additional information can be found at www.cognivue.com. Cognivue Support personnel can be reached from 9-5pm EST at 1-585-433-2992 or support@cognivue.com.

UM-450-G Cognivue Portal User Manual, 2/14/2023



7911 Rae Blvd | Victor, NY 14564 | 585.203.1969

1.0 About the Cognivue Portal

Cognivue Clarity® Device is a computerized cognitive assessment aid that objectively, quantitatively and reliably evaluate cognitive function. Cognivue Thrive® Device is a computerized cognitive assessment aid that objectively, quantitatively and reliably screens out cognitive impairment. After a user completes a series of subtests, both devices produce easy-to-interpret reports with scores of various cognitive domains. The test results and associated reports are stored on the Cognivue device used to take the test.

The Cognivue Portal enables you to access your Cognivue Clarity Device or Cognivue Thrive Device data from anywhere you have access to the internet. In addition, the Portal makes it easier for you to search for and find reports. Specifically the Portal allows authorized users in your organization to:


- Search for patients who have taken a Cognivue test on your devices
- View a summary of the test results available for that patient
- View reports for that patient
- View a trend of the scores for that patient
- List and view all reports across all your patients and devices
- Download and share reports

The Portal is to be used by subscribed Cognivue Customers only and is not intended to provide direct patient access to Cognivue Test Reports. Its purpose is solely to transfer, store and display results. It is the responsibility of the Customer to provide patient access to Patient Data in accordance with their data privacy and security procedures and guidelines. This may include procedures such as giving a printed copy to a patient or sending secure email (outside of the Cognivue product).

For information about the Cognivue Clarity Device and Thrive Device, including a description of the intended use of the devices and their test results, please refer to their respective User Manual.

2.0 Portal Description

The Cognivue Portal uses secure cloud-based storage to enable access across all your devices (see the figure below).



As shown, you can use any modern web browser to access the portal. Access to the portal is controlled separately from access on the device. In order to access data via the portal the following must occur:

- 1) The device must be configured to sync with the portal. Each device needs to be configured separately.
- 2) Cognivue Device users must separately be given access to the portal.
- 3) Portal users need to be assigned to devices that they have access to.

Cognivue implements a management approved and funded Information Security and Privacy Program that is HIPAA Compliant and ISO 27001 aligned. This program has been implemented to help ensure the confidentiality, integrity, availability, and privacy of its clients and partners. The program includes administrative, physical and technical controls that include, but are not limited to, the following:

1. Risk Management and Assessments
2. Data Management and Classification
3. Asset and Vendor Management, including BAAs in place for all critical vendors
4. Physical Security Controls
5. Workforce Security Controls
6. Identity and access Management Controls
7. Software Development Lifecycle and Change Management
8. Infrastructure Security Controls
9. Contingency Planning, Disaster Recovery and Business Continuity
10. Privacy and Data Use Controls

3.0 Recommendations for Use

To provide an efficient, secure use of the portal, Cognivue recommends the following practices.

- 1) Follow your organizational procedures for protecting patient security, privacy, and access to reports.
- 2) The Portal is developed for customer staff only. Patients should not have direct access or a login to the portal.
- 3) Implement Minimum Necessary access to accounts. Device Administrators do not need to be Portal Administrators. It is not recommended to give Portal Administrator access to a user unless they are intended to manage your portal accounts and settings.
- 4) Cognivue reports should be shared with patients using HIPAA compliant and secure procedures maintained by the customer organization. This may include secure email, paper printout or EMR integrations.
- 5) Passwords for the portal should be long and complex and be developed according to the customer organization's security policies. Where possible they should be random and stored in a password manager.


4.0 Getting Started & User Management

4.1 Configuring a Cognivue Device to Sync to the Portal


An Administrator can enable this functionality by going to the Portal menu in the main window and clicking on settings.

To enable the device to sync to the portal, the Device Administrator first selects “Portal” and “Settings” from the main window.


Main Device Window – Device Administration Menu – Clarity




Main Device Window – Device Administration Menu - Thrive



In the portal settings window, click on the “Enable Sync” button(circled below).



Fill in the information and click the save button. This will create the initial administrator account in the portal.





The screenshot shows the "Portal Settings" window with the same red message box. Below the message box is an "Enable Sync" button. Below that is a text prompt: "You must create a portal account in order to use the portal. Please fill out the fields below, then click 'Save'." There are four required text input fields: "*Email Address", "*Confirm Email", "*First Name", and "*Last Name". At the bottom right, there is a "Save" button circled in red. A legend at the bottom left indicates that the asterisk (*) denotes a required field.

By default, the data sync will begin upon completion of the next patient test. If an administrator wants to force the sync to occur immediately, they can select the “Force Sync” option from the


Portal menu to send data from the device to the Portal system. At the main window, click the “Portal” menu and select “Force Sync”.

Force Sync – Clarity



Force Sync - Thrive

A pop-up window will allow you to confirm or cancel the start of the data sync.



The next step is to setup the Portal Accounts.


4.2 Setting up the First Portal Account

The initial portal administrator account is created with the information provided by the device administrator, when enabling the sync in the portal settings window.

4.3 Granting Device Users to Access the Portal (Provisioning Users)

Device Users don't automatically have access to the Portal - instead, they must be granted access to the Portal by a process called "provisioning". To do this, a Portal Administrator uses the **Provision Users** menu on the Portal to grant device users access to the portal. The process for provisioning users is provided below.

1. Log into the [Portal](#) as an administrator.
2. Navigate to the Admin drop-down button, then select Provision Users:



Portal Provision Users Window

3. Find the user you would like to give access to, then click Provision User:


The screenshot shows the 'Provision User' form. It has a title bar with 'Provision User' and a close button. The form contains the following fields: Email (operator2@gmail.com), First Name (Operator), Last Name (Two), Portal Role (User), and a section for selecting devices with a 'Select All' button. At the bottom, there are 'Close' and 'Submit' buttons.

4. Ensure the user's information is correct
5. Assign the correct user role
6. The user will have access to information from the devices selected from their account. Users with the Administrator role will have access to information from all testing devices for their account. Administrators also have the ability to manage and provision Users.
7. When finished, click the Submit button
8. The user will receive an email asking them to finish setting up their Portal account. The email will contain a link to the Portal and a temporary password.

4.4 Adding a New User that is not a Device User to the Portal

A user can be created directly in the Portal, rather than provisioned from a user on the device. This user will be associated with the same account as the admin who created it.

1. From the Portal Admin Page – Manage Users click the Create User button.



2. Fill in the information for the User:

The 'Create User' form is displayed. It includes the following fields and options:

- Email:
- First Name:
- Last Name:
- Portal Role:
- Select which devices this user can access:
 -
 -
 -

At the bottom right, there are 'Close' and 'Submit' buttons.

Create User Window

3. Assign the desired user role for the user. It is not recommended to give Administrator access to a user unless they are intended to manage your portal accounts and settings.

4. The user will have access to information from the devices selected from their account. Users with the Administrator role will have access to information from all testing devices for the account. Devices deactivated and then reactivated, will be selectable if the device has been used to run tests for your Account, previously.
5. When finished, click the Submit button.
6. The user will receive an email asking them to finish setting up their account.

4.5 Editing a user

The characteristics of existing users can be modified if the 'Edit User' button is selected on the Portal Admin page.

The screenshot shows the Cognivue Portal Admin interface. At the top, there are navigation buttons for Home, Reports, Export Data, and Admin. A dropdown menu for Admin is open, showing options like Manage Users, Provision Users, Share History, and Manage Devices. The main content area is titled 'Manage Existing Users' and includes a search bar and a 'Create User' button. Below this is a table with the following data:

First Name	Last Name	Email	Portal Role	Confirmed Email	Two-Factor Enabled	Active	
Portal	Admin	cogportaldev+pa01@gmail.com	Administrator	✓		✓	Edit User
First	Admin	cogportaldev+firstadmin@gmail.com	Administrator	✓			Edit User
First	Operator	cogportaldev+fo01@gmail.com	User	✓		✓	Edit User

The 'Edit User' button for the user 'First Operator' is circled in red. Below the table, there is a note: 'Table containing list of portal user accounts and information about them'. At the bottom right, there are 'Previous', '1', and 'Next' navigation controls.

The window shown below can then be used to modify the appropriate user data. The button on the bottom of the window can be used to either deactivate or activate users.

The 'Edit User' window is a modal form for editing user information. It contains the following fields and controls:


- Email: cogportaldev+fo01@gmail.com
- First Name: First
- Last Name: Operator
- Portal Role: User (dropdown menu)
- Select which devices this user can access: Select All button and a list of selected devices (COG1CAVT2 \Portal Device).
- Deactivate User button (blue)
- Close and Save buttons (bottom right)

Edit User Window

5.0 Portal Functionality

5.1 Searching for Patients & The Patient Dashboard

When first logging into the site, or after clicking the New Search button from the Patient Information Dashboard, the patient search box will be available:




The screenshot shows a 'Patient Search' window with the following elements:

- Cognivue logo
- Section title: Patient Search
- Input fields: First Name, Middle Initial, Last Name *, Date of Birth *
- Date of Birth field includes: mm/dd/yyyy, a calendar icon, and a Clear button
- Sex selection: Male, Female, Not Specified (radio buttons)
- Search button


Patient Search Window

1. The minimum information required to search is either Date of Birth or Last Name, though providing more information will help make finding the right person easier.
2. If your search produces multiple results, you will be asked to select a single person from a list:



Select Patient Window

- Once you've selected a patient, you will be taken to the patient dashboard that contains an overview of that patient's most recent test results.




The left panels of the dashboard provides patient information, the ability to start a new search, and some useful links for contacting Cognivue.

In the center panel, the latest test score and test score trend are displayed.

In the right panel, a listing of test scores and recent reports are provided. If you select the links on this panel, you will be able to share or download the selected report.

5.2 The Reports Page

When the Reports Button at the top of the window is selected, a listing of all reports on all the devices that you have access to will be provided. This listing can be sorted by any of the columns as well as searched by patient name or test report range.




The screenshot displays the Cognivue portal interface. At the top, there is a navigation bar with buttons for Home, Reports, Export Data, and Admin. A user profile for 'Hello Portal User!' is visible on the right. Below the navigation bar, there is a search bar and a 'Test Date Range' input field. The main content area features a table with 11 rows of report data. Each row includes a checkbox for selection, a name, date of birth, sex, test date, and two action icons (a cloud and a page with an arrow). Below the table, there is a 'Download Selected' button and a pagination control showing 'Previous 1 2 3 Next'.

<input type="checkbox"/>	Name	Date of Birth	Sex	Test Date	Actions
<input type="checkbox"/>	Some, Name X	December 6, 1994	Male	April 17, 2022 8:08 PM	
<input type="checkbox"/>	Castro, Honey A	January 9, 1988	Female	August 6, 2021 6:21 PM	
<input type="checkbox"/>	Castro, Lily Y	January 8, 1987	Female	July 5, 2021 6:21 PM	
<input type="checkbox"/>	Chan, Sunny N	January 7, 1986	Male	March 8, 2021 6:21 PM	
<input type="checkbox"/>	Castro, Simon S	January 6, 1985	Male	March 7, 2021 6:21 PM	
<input type="checkbox"/>	Castro, Daisy D	January 5, 1984	Female	March 6, 2021 6:21 PM	
<input type="checkbox"/>	Chan, Diego D	January 4, 1983	Male	March 5, 2021 6:21 PM	
<input type="checkbox"/>	Chan, Luna L	January 3, 1982	Male	February 5, 2021 6:21 PM	
<input type="checkbox"/>	Chan, William A	January 2, 1981	Male	January 5, 2021 6:21 PM	
<input type="checkbox"/>	Chan, Sunny D.	January 1, 1980	Female	January 5, 2021 6:11 PM	

Individual reports can be downloaded by clicking on the cloud icon in the right-most column. Reports may also be emailed to recipients with the page-with-arrow icon in the right-most column. You can also select multiple reports in the first column and download them at the same time using the “Download Selected” button at the bottom of the screen.


5.3 Downloading Test Data

An Administrator Account can export all of their patient test data. Click on the Export Data option on the navigation bar, to go to the Export Data page.



Home Page

On the Export Data page, click on the 'Export to CSV' button. This will download all patient test results for the administrator's account, into a .zip file containing the patient test results in a .csv file. If there are test results for both Clarity and Thrive devices on the Administrator's account, then there will be one .csv file for any Clarity test results, and one .csv file for any Thrive test results.




Export Data Page

5.4 History of Shared Reports

An Administrator Account is able to view the history of shared reports for their patients and devices. To navigate there, click on the Admin drop-down menu and select the Share History option.



Admin Menu – Share History



Share History Page

The Administrator Account is able to navigate to the patient dashboard by clicking on the patient’s name.


As with the report page, the Account may select to show 10, 25, 50, or 100 records per page, and may navigate between the pages of the table by using the buttons below the table.

The ‘Downloaded Date’ column reflects the state of the shared report link. If a Portal Account has shared a report multiple times, there will be one entry on this table for each time that report has been shared. If the Downloaded Date is empty, then the shared report has not been downloaded and the link is still valid for use. If ‘Link Expired’ appears, then the shared report was not downloaded and the link has expired. When the shared report is downloaded a timestamp will appear in this field.

6.0 Devices Page

6.1 Accessing the Devices Page


The User is able to review their devices by accessing the devices page. The User may access the Page by clicking on the Admin button and selecting 'Manage Devices' from the list. This will bring the User to the Devices Page.



Manage Devices Button

6.1.1 Adding a Nickname to a Device

The User can add or change the portal nickname for a device by selecting the Edit button in the Actions column, for the device. This will open a window with the device information and a field to add or change the nickname.



Edit Button

Edit Device ✕

Serial Number:

Device Nickname:

Is Device Active?

Activation Date:

Device Deactivation Date:

Close Save


Edit Device Modal

Fill in the Device Nickname field and click the Save button, to add or change the portal nickname assigned to the device.


7.0 Account Management

7.1 Accessing the Account Management Page

The User is able to manage various aspects of their account by accessing the account management page. The User may access the page by clicking on the button in the top right of any page, which says 'Hello' followed by their username, and select 'Account. This will bring the User to the Account Management page.



Account Management Button



Account Management Page

7.2 Profile Page

On this page the User is able to change their first name and/or last name. To do so, edit the contents of either field and click the Save button.

Manage your account

Profile

Username
cogportaldev+portaladmin@gmail.com

First name
Portal

Last name
Admin

Save

Profile Page

7.3 Password Page

On this page the User is able to change their password for logging into the portal. To do so, enter the current password in the first field, and the new password in the second and third fields. Passwords must contain at least 1 uppercase letter, 1 lowercase letter, 1 number, and 1 special character, and must be a minimum of 8 characters.

Manage your account

Change password

Current password

New password


Confirm new password

Update password

Password Page

7.4 Personal Data Page


On this page the User is able to delete their portal user account and associated data. The deleted data includes the phone number and email associated with the User account.



Personal Data Page

7.5 Two-Factor Authentication

In order to provide better security for Users, a User is able to set up two-factor authentication with an authenticator app and/or a phone number where they can receive verification codes through texts.



Two-Factor Authentication Page with no 2FA set up

Manage your account



The screenshot shows the 'Manage your account' page with a sidebar on the left containing links for Profile, Password, Personal data, and Two-Factor Authentication (which is highlighted). A green notification banner at the top states 'Your authenticator app has been verified.' The main content area is titled 'Two-Factor Authentication' and includes the text 'Two-factor authentication is enabled for your Portal account.' Below this are buttons for 'Reset Recovery Codes', 'Add Authenticator App', 'Reset Authenticator App', and 'Disable Authenticator App'. The 'Phone Number' section indicates 'SMS authentication is enabled.' and includes buttons for 'Change Phone Number' and 'Disable SMS Authentication'. At the bottom, there is a 'Disable 2FA' button.

Two-Factor Authentication with multiple methods configured

7.5.1 Setup Authenticator App

If you select the Setup Authenticator App button, you can connect an authenticator app that you have, to your user account. Follow the instructions of the authenticator app that you have installed, and the on-screen instructions on this page, to connect your authenticator app with your user account.

Manage your account

The screenshot shows the 'Configure authenticator app' page. It includes a sidebar with links for Profile, Password, Personal data, and Two-Factor Authentication (highlighted). The main content area is titled 'Configure authenticator app' and provides instructions: 'To use an authenticator app go through the following steps:'. Step 1: 'Download a two-factor authenticator app like Microsoft Authenticator for Windows Phone, Android and iOS or Google Authenticator for Android and iOS.' Step 2: 'Scan the QR Code or enter this key:   authenticator app. Spaces and casing do not matter.' Step 3: 'Once you have scanned the QR code or input the key above, your authenticator app will provide you with a unique code. Enter the code in the confirmation box below.' Below the instructions is a 'Verification Code' input field and a 'Verify' button.

Setup Authenticator App

After connecting to your authenticator app, enter the verification code into the field and click the 'Verify' button, to complete the setup for the User account.

7.5.2 Setup SMS (texting)

If the User selects the Setup SMS button, they can configure a phone number to receive codes via SMS (Texting.) Follow the on-screen instructions to complete the configuration. Only one phone number may be configured for SMS for the User account.

Manage your account

Setup Authentication through SMS

Enter the phone number that you will use to receive authentication codes via SMS. After saving, you will receive a text message with a verification code. You will then be prompted to enter the code, which must be done before you can login with SMS authentication. If you do not have another form of two-factor authentication set up, Two Factor Authentication will not be enabled until you enter the verification code.

Phone Number

Confirm

Setup SMS Page

Manage your account

Verify Your Phone Number

A verification code has been sent to [redacted]. Enter [redacted] now.

SMS authentication will not be enabled until you verify your phone number. If you leave this page and do not enter the verification code you received, SMS authentication will not be enabled, and you will have to go through the setup process again to enable it.

Verification Code

Verify Resend Code

Verify Phone Number Page

7.5.3 Recovery Codes

Recovery codes are single-use codes that can be used to verify the User account. These are intended to be used in the event that the User does not have access to the authentication method(s) that have been configured for their User account. For example, the User does not have their mobile device with the authenticator app, available.

When the User initially configures either the authenticator app or SMS, for two-factor authentication, 10 recovery codes will be automatically generated and displayed to the screen. The User will be unable to retrieve these codes once they leave that page.

The User can manually generate recovery codes by clicking on the Reset Recovery Codes button. This will take the User to the Generate Recovery Codes page. Follow the on-screen instructions to generate 10 recovery codes and display them to the screen. The User will be unable to retrieve them once they leave the page.

Manage your account

The screenshot shows the 'Manage your account' page with a navigation menu on the left containing 'Profile', 'Password', 'Personal data', and 'Two-Factor Authentication'. A green notification bar at the top states 'Your authenticator app has been verified.' The main content area is titled 'Two-Factor Authentication' and includes the text 'Two-factor authentication is enabled for your Portal account.' Below this, the 'Reset Recovery Codes' button is circled in red. Further down, the 'Authenticator App' section shows 'Authenticator app is enabled.' with buttons for 'Add Authenticator App', 'Reset Authenticator App', and 'Disable Authenticator App'. The 'Phone Number' section shows 'SMS authentication is enabled.' with buttons for 'Change Phone Number' and 'Disable SMS Authentication'. At the bottom, there is a 'Disable ZFA' button.

Two-factor Authentication Page; Reset Recovery Codes button

Manage your account

The screenshot shows the 'Manage your account' page with a navigation menu on the left containing 'Profile', 'Password', 'Personal data', and 'Two-Factor Authentication'. The main content area is titled 'Generate two-factor authentication recovery codes'. A yellow warning box contains the text: 'Put these codes in a safe place. If you lose your device and don't have the recovery codes you will lose access to your account. Generating new recovery codes does not change the keys used in authenticator apps. If you wish to change the key used in an authenticator app you should reset your authenticator keys.' Below the warning box is a 'Generate Recovery Codes' button.

Generate Recovery Codes Page

Manage your account

The screenshot shows the 'Manage your account' page with a navigation menu on the left containing 'Profile', 'Password', 'Personal data', and 'Two-Factor Authentication'. A green notification bar at the top states 'You have generated new recovery codes.' The main content area is titled 'Recovery codes'. A yellow warning box contains the text: 'Put these codes in a safe place. If you lose your device and don't have the recovery codes you will lose access to your account.' Below the warning box, five recovery codes are listed: 542ab8ce e3d258bd, 73bb81ad 276c6185, e9f9d1be 59aa8a97, f71512a2 06a2879b, and 7a914aa2 69ad4d1c.

Recovery Codes Page

7.5.4 Resetting Authenticator App

The User can reset the authenticator app to disconnect it and allow them to connect to a different authenticator app. To reset the authenticator app security on their User account, the User must click on the Reset Authenticator App button.

Manage your account

Profile
Password
Personal data
Two-Factor Authentication

Your authenticator app has been verified. ×

Two-Factor Authentication

Two-factor authentication is enabled for your Portal account.

[Reset Recovery Codes](#)

Authenticator App

Authenticator app is enabled.

[Add Authenticator App](#) [Reset Authenticator App](#) [Disable Authenticator App](#)

Phone Number

SMS authentication is enabled.

[Change Phone Number](#) [Disable SMS Authentication](#)

[Disable 2FA](#)

Two-Factor Authentication Page; Reset Authenticator App button

Follow the instructions on the Reset Authenticator App page, in order to reset the authenticator app for the User account.

Manage your account

Profile
Password
Personal data
Two-Factor Authentication

Reset authenticator key

If you reset your authenticator key you will not be able to sign in with your authenticator app until you reconfigure it.

You will be redirected to the authenticator app setup page where you can reconfigure your authenticator app with the new key. If you do not have SMS authentication set up, this action will disable two-factor authentication until you verify your authenticator app again, or set up SMS authentication.

[Reset authenticator key](#)

Reset Authenticator App Page

Manage your account


Profile
Password
Personal data
Two-Factor Authentication

Your authenticator app key has been reset, you will need to configure your authenticator app using the new key. You will not be able to log in with your authenticator app until you do. If you do not have SMS authentication enabled, two-factor authentication will be disabled until you reconfigure your authenticator app or set up SMS authentication. ×

Configure authenticator app

To use an authenticator app go through the following steps:

1. Download a two-factor authenticator app like Microsoft Authenticator for [Windows Phone, Android](#) and [iOS](#) or Google Authenticator for [Android](#) and [iOS](#).
2. Scan the QR Code or enter this key `455F-3466-3223` authenticator app. Spaces and casing do not matter.



3. Once you have scanned the QR code or input the key above, your authenticator app will provide you with a unique code. Enter the code in the confirmation box below.

Verification Code

[Verify](#)


Setup Authenticator App Page after resetting the authenticator app

7.5.5 Disabling Two-Factor Authentication

The User may disable two-factor authentication on their account, either for a specific method, or completely. It is recommended that the User maintain at least one configured method for two-factor authentication.


7.5.5.1 Disabling Authenticator App 2FA

To disable two-factor authentication for the authenticator app only, click on the Disable Authenticator App Authentication button.



Two-Factor Authentication Page; Disable Authenticator App button


Follow the on-screen instructions to disable the authenticator app for the User.



Disable Authenticator App Page


7.5.5.2 Disable SMS 2FA

To disable two-factor authentication for SMS only, click on the Disable SMS Authentication button.



Two-Factor Authentication Page; Disable SMS Authentication button

Follow the on-screen instructions to disable SMS authentication for the User.



Disable SMS Authentication Page

7.5.5.3 Disable Two-Factor Authentication

To disable all two-factor authentication for their user account, click the Disable 2FA Button. **It is recommended that a User have at least one method of authentication configured for their account.**

Manage your account

The screenshot shows the 'Manage your account' page with a navigation menu on the left containing 'Profile', 'Password', 'Personal data', and 'Two-Factor Authentication'. A green notification bar at the top states 'Your authenticator app has been verified.' The main content area is titled 'Two-Factor Authentication' and includes the text 'Two-factor authentication is enabled for your Portal account.' Below this are buttons for 'Reset Recovery Codes', 'Add Authenticator App', 'Reset Authenticator App', and 'Disable Authenticator App'. The 'Authenticator App' section shows 'Authenticator app is enabled.' and buttons for 'Add Authenticator App', 'Reset Authenticator App', and 'Disable Authenticator App'. The 'Phone Number' section shows 'SMS authentication is enabled.' and buttons for 'Change Phone Number' and 'Disable SMS Authentication'. A red circle highlights the 'Disable 2FA' button at the bottom of the page.

Two-Factor Authentication Page; Disable 2FA button

Follow the on-screen instructions to disable two-factor authentication for all methods, on the Useraccount.

Manage your account


The screenshot shows the 'Manage your account' page with a navigation menu on the left containing 'Profile', 'Password', 'Personal data', and 'Two-Factor Authentication'. The main content area is titled 'Disable Two-Factor Authentication' and includes a yellow warning box with the text 'This action disables two-factor authentication.' and 'Disabling two-factor authentication will reset your authenticator key and phone number. In order to re-enable two-factor authentication, you will have to go through the set-up process again for either an authenticator app or SMS authentication. It is recommended that you protect your account with two-factor authentication. Disabling it may leave your account more vulnerable.' Below the warning box is a red button labeled 'Disable two-factor authentication'.

Disable 2FA Page

8.0 Portal Support

8.1 Accessing the Portal Support Page

The User is able to submit requests related to the Portal to Cognivue. The User can submit requests for portal support, and submit requests for enhancements to the portal. The User may access the page by clicking on the button in the top right of any page, which says 'Hello' followed by their username, and select 'Support'. This will bring the User to the Portal Support page.



Support Button

Cognivue® Portal User Manual

The screenshot shows a web browser window with the URL `dev-portal.cognivue.com/SupportRequest/RequestSupport`. The page title is "Support Request" and the subtitle is "Submit Support Request". The form includes the following fields and options:

- Type of Request:** Radio buttons for "Support" and "Enhancement".
- What product is your request for?:** A dropdown menu with "Portal" selected.
- Account Name:** A text input field containing "Aaron Testing Account".
- Name of Requester:** An empty text input field.
- Email:** A text input field containing "cogportaldev+pa01@gmail.com".
- Phone Number (Optional):** An empty text input field.
- How can we help you today?:** A large text area for describing the request.
- May we reach out for further details?:** A checked checkbox.
- Would you like us to give you status updates?:** A checked checkbox.
- Submit:** A button at the bottom of the form.

Support Request Form

The screenshot shows the same web browser window as above, but with the "Type of Request" set to "Enhancement". The form includes the following fields and options:

- Type of Request:** Radio buttons for "Support" and "Enhancement".
- What product is your request for?:** A dropdown menu with "Device" selected.
- Serial number:** A dropdown menu with "Select a Serial Number" and a "Search SN" input field below it.
- Account Name:** A text input field containing "Aaron Testing Account".
- Name of Requester:** An empty text input field.
- Email:** A text input field containing "cogportaldev+pa01@gmail.com".
- Phone Number (Optional):** An empty text input field.
- How can we help you today?:** A large text area for describing the request.
- May we reach out for further details?:** A checked checkbox.
- Would you like us to give you status updates?:** A checked checkbox.
- Submit:** A button at the bottom of the form.

Support Request Form – Device Requests

8.2 Submit a Request

The User can toggle between Support and Enhancement requests by clicking on the buttons labeled 'Support' and 'Enhancement' on the Support page.

8.2.1 Support Request

The User can submit requests for portal support on the Submit Support Request page.

- *What product is your request for?* The User making the request should select from the dropdown menu whether their request is for the portal, a device, or some other unclassified Cognivue service or product.
- *Serial Number (for Device Support Requests Only):* The User selects the serial number of the device that they are submitting the request for. This field only appears when the User selects the Device option on the Support Page.
- *Search SN:* The User may type into this field to search for a specific device serialnumber. The serial number dropdown menu will try and select the serialnumber that best matches the current value of the field, as the User types.
- *Account Name:* The name of the User's account that their devices are classified under.
- *Name of Requester:* The User making the request should enter their name here
- *Email:* The User making the request should enter the email address where they would like any correspondence sent to; a copy of the support request will be sent to this email address. This field will be initially populated with the email of the Portal User that is logged in.
- *Phone Number:* The User may optionally include a phone number that they may be contacted at.
- *How can we help you today?:* This is where the User should enter what portal issue they are experiencing including details on what happened, what was going on at the time, and any actions performed.
- *May we reach out for further details?:* The User should check this box if they would like Cognivue Support to reach out to them for further information and clarification.
- *Would you like us to give you status updates?:* The User should check this box if they would like Cognivue Support to provide them with updates on the status of their Support Request.

Cognivue® Portal User Manual

The screenshot shows the 'Support Request' form on the Cognivue portal. The page title is 'Support Request' and the subtitle is 'Submit Enhancement Request'. The 'Type of Request' section has two buttons: 'Support' (highlighted in blue) and 'Enhancement'. Below this is a dropdown menu for 'What product is your request for?' with 'Portal' selected. The form also includes fields for 'Account Name' (Aaron Testing Account), 'Name of Requester', 'Email' (cogportaldev+pa01@gmail.com), and 'Phone Number (Optional)'. There are two text areas for 'How can we help you today?' and 'What is the intended use of the Request?'.

This screenshot shows the bottom portion of the 'Support Request' form. It includes the 'How can we help you today?' and 'What is the intended use of the Request?' text areas. Below these are two checkboxes: 'May we reach out for further details?' and 'Would you like us to give you status updates?', both of which are checked. At the bottom of the form is a 'Submit' button.

Enhancement Request Form

8.2.2 Enhancement Request

The User can submit requests for portal enhancements on the Submit Enhancement Request page.

- *What product is your request for?* The User making the request should select from the dropdown menu whether their request is for the portal, a device, or some other unclassified Cognivue service or product.
- *Account Name:* The name of the User's account that their devices are classified under.
- *Name of Requester:* The User making the request should enter their name here
- *Email:* The User making the request should enter the email address where they would like any correspondence sent to; a copy of the support request will be sent to this email address. This field will be initially populated with the email of the Portal User that is logged in.
- *Phone Number:* The User may optionally include a phone number that they may be contacted at.
- *How can we help you today?:* This is where the User should enter what they would like to see changed or added in the portal.
- *What is the intended use of the Request?:* This is where the User should describe what the request is intended for. The User should provide examples and scenarios for application.
- *May we reach out for further details?:* The User should check this box if they would like Cognivue Support to reach out to them for further information and clarification.
- *Would you like us to give you status updates?:* The User should check this box if they would like Cognivue Support to provide them with updates on the status of their Support Request.